




# Integrating Four Core Technologies TO SHAPE THE FUTURE OF MEDICINE

ASAHI® IS THE INDUSTRY LEADER IN GUIDE WIRE TECHNOLOGY. INTEGRATION OF 4 CORE PROPRIETARY TECHNOLOGIES PROVIDES PURPOSEFUL OPTIMIZATION OF GUIDE WIRE PERFORMANCE FOR BOTH ROUTINE AND COMPLEX PROCEDURES.



\*Available in select ASAHI guide wires

## ORDERING INFORMATION

PERIPHERAL GUIDE WIRES					
PRODUCT	CATALOG NO.	DIAMETER	LENGTH	ACT ONE	TIP SHAPE
Regalia XS 1.0 0.014 inch	PAGP140000	0.014 inch (0.36mm)	180cm	—	Straight
	PAGP140300	0.014 inch (0.36mm)	300cm	—	Straight
ASAHI Gladius 0.014 inch	PPW14R100S	0.014 inch (0.36mm)	200cm	ACT ONE	Straight
	PPW14R300S	0.014 inch (0.36mm)	300cm	ACT ONE	Straight
ASAHI Halberd 0.014 inch	PHW14R101S	0.014 inch (0.36mm)	200cm	ACT ONE	Straight
	PHW14R301S	0.014 inch (0.36mm)	300cm	ACT ONE	Straight
Astato XS 20 0.014 inch	PAGH143092	0.014 inch (0.36mm)	180cm	—	Straight
	PAGH143392	0.014 inch (0.36mm)	300cm	—	Straight
Astato XS 40 0.014 inch	PAGHW143094	0.014 inch (0.36mm)	200cm	—	Straight
	PAGHW143394	0.014 inch (0.36mm)	300cm	—	Straight
ASAHI Gladius 0.018 inch	PPW18R100S	0.018 inch (0.45mm)	200cm	ACT ONE	Straight
	PPW18R300S	0.018 inch (0.45mm)	300cm	ACT ONE	Straight
Treasure Floppy 0.018 inch	PAGH18M072	0.018 inch (0.45mm)	190cm	—	Straight
	PAGH18M372	0.018 inch (0.45mm)	300cm	—	Straight
ASAHI Gaia PV 0.018 inch	PHW18R102P	0.018 inch (0.45mm)	200cm	ACT ONE	Pre-shaped
	PHW18R302P	0.018 inch (0.45mm)	300cm	ACT ONE	Pre-shaped
Treasure 12 0.018 inch	PAGH18M070	0.018 inch (0.45mm)	180cm	—	Straight
	PAGH18M370	0.018 inch (0.45mm)	300cm	—	Straight
ASAHI Halberd 0.018 inch	PHW18R101S	0.018 inch (0.45mm)	200cm	ACT ONE	Straight
	PHW18R301S	0.018 inch (0.45mm)	300cm	ACT ONE	Straight
Astato 30 0.018 inch	PAGH18M071	0.018 inch (0.45mm)	180cm	—	Straight
	PAGH18M371	0.018 inch (0.45mm)	300cm	—	Straight

**Indications for Use** The ASAHI INTECC peripheral guide wires are intended to facilitate the placement and exchange of diagnostic and therapeutic devices during intravascular procedures. These devices are intended for peripheral vascular use only.

©2018 ASAHI INTECC CO., LTD.

“ASAHI”, “ACT ONE”, “Regalia”, “ASAHI Gladius” and “Halberd” are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Some materials / graphics were provided by Vascular Perspectives – [www.vascularperspectives.com](http://www.vascularperspectives.com).



ASAHI INTECC USA, INC.  
Toll Free: 855-286-9473 Fax: 949-756-8146  
WWW.ASAHI-INTECC-US.COM

ASAHI Peripheral Guide Wires/ver.2/L/AMC-K18231  
MKT AIU 202 Rev. A




# PERIPHERAL GUIDE WIRES

MASTER THE EVERYDAY

ASAHI INTECC USA, INC.

### 0.014 Inch Guide Wires



### BELOW THE KNEE

#### 0.014 inch (0.36mm)


Product	Tube Clip	TIP LOAD
<b>Regalia<sup>®</sup> XS 1.0</b> Sub-total / Tortuous		1gf
<b>ASAHI Gladius<sup>®</sup></b> 0.014 inch Workhorse Plus		3gf
<b>ASAHI Halberd<sup>®</sup></b> 0.014 inch Controlled Navigation		12gf
<b>Aстато<sup>®</sup> XS 20</b> 0.014 inch Highly Calcified Lesion		20gf
<b>Aстато XS 40</b> 0.014 inch Highly Calcified Lesion		40gf

#### 0.018 inch (0.45mm)

<b>ASAHI Gladius</b> 0.018 inch Workhorse		4gf
<b>Treasure<sup>®</sup> 12</b> Controlled Drilling		12gf




### WIRE COATINGS LEGEND



\*Coated with SLIP-COAT<sup>®</sup> coating

### 0.018 Inch Guide Wires



### ABOVE THE KNEE

#### 0.018 inch (0.45mm)

Product	Tube Clip	TIP LOAD
<b>ASAHI Gladius</b> 0.018 inch Workhorse		4gf
<b>Treasure Floppy<sup>®</sup></b> 0.018 inch Workhorse		4gf
<b>Treasure 12</b> Controlled Drilling		12gf
<b>ASAHI Gaia<sup>®</sup> PV</b> 0.018 inch Complex Lesion		7.5gf
<b>ASAHI Halberd</b> 0.018 inch Controlled Navigation		12gf
<b>Aстато 30</b> 0.018 inch Highly Calcified Lesion		30gf

#### 0.014 inch (0.36mm)

<b>Aстато XS 20</b> 0.014 inch Highly Calcified Lesion		20gf
<b>Aстато XS 40</b> 0.014 inch Highly Calcified Lesion		40gf