

Captivator™ II

Single-Use Snares



Captivator™ II

Single-Use Snares

Captivator II Snares are the first line of stiff and rounded snares, available in multiple sizes with both a hot and cold snaring indication.

Features and Benefits

- | | |
|------------------------------|--|
| Rounded Shape | <ul style="list-style-type: none">• Short distance from scope tip is designed for tip control / visualization• Better capture of flat lesions |
| Stiff Braided Wire | <ul style="list-style-type: none">• Designed to provide grip / traction on flat polyps, Control of wire |
| Hot / Cold Indication | <ul style="list-style-type: none">• Allows clinical flexibility |
| Smaller Loop Options | <ul style="list-style-type: none">• Diminutive polyps• Residual tissue during EMR |

PRODUCT SPECIFICATIONS

Loop Diameters

- 10mm
- 15mm
- 20mm
- 25mm
- 33mm

Sheath Diameter

- 2.4mm

Working Length

- 240cm

Working Channel

- 2.8mm



ORDERING INFORMATION

Captivator II Single-Use Snares

| Order Number | Description | Loop Width (mm) | Sheath OD (mm) | Working Length (cm) | Min. Working Channel (mm) | Packaging |
|--------------|-----------------------------|-----------------|----------------|---------------------|---------------------------|-----------|
| M00561221 | 10mm Rounded – Stiff | 10mm | 2.4 | 240 | 2.8 | Box 10 |
| M00561222 | 10mm Rounded – Stiff | 10mm | 2.4 | 240 | 2.8 | Box 20 |
| M00561231 | 15mm Rounded – Stiff | 15mm | 2.4 | 240 | 2.8 | Box 10 |
| M00561232 | 15mm Rounded – Stiff | 15mm | 2.4 | 240 | 2.8 | Box 20 |
| M00561241 | 20mm Rounded – Stiff | 20mm | 2.4 | 240 | 2.8 | Box 10 |
| M00561242 | 20mm Rounded – Stiff | 20mm | 2.4 | 240 | 2.8 | Box 20 |
| M00561191 | 25mm Rounded – Stiff | 25mm | 2.4 | 240 | 2.8 | Box 10 |
| M00561192 | 25mm Rounded – Stiff | 25mm | 2.4 | 240 | 2.8 | Box 20 |
| M00561291 | Extra Large Rounded – Stiff | 33mm | 2.4 | 240 | 2.8 | Box 10 |

Boston Scientific

Advancing science for life™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/endo-resources

Ordering Information
1.888.272.1001

©2014 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-214401-AA May 2014

Captivator is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.